

## **Kidevo** mini and mini.t

Rigid frame wheelchairs for kids.



Feature	Kidevo mini	Kidevo mini.t
Max. user weight	50 kg	50 kg
Seat width range	16 – 32 cm	16 – 32 cm
Seat depth range	16 – 35 cm	16 – 35 cm
Tilt	–	-5 to + 40°
Drive wheel sizes	20, 22 and 24"	20, 22 and 24"
Crash tested (acc. ISO 7176-19)	✓	✓



### Easy to use brakes

- All *Kidevo* wheelchairs can be chosen with different friendly brakes such as the drum brake. Combinations of brakes are also possible.



### Functional seat tilt

- The *Kidevo mini.t* has a functional and space-saving seat tilt without changing the center of gravity and thereby, maintaining an optimal driving behavior.




### Robust foot plate

- For an optimal foot positioning, the *Kidevo* wheelchairs come with a robust and flat one piece foot plate. In addition, there is a high variety of more positioning options to choose from.



### High variety of rear wheels

- Each model of the *Kidevo* provides three different rear wheel sizes in combination with various types of tires and handrims.



### Custom-made seat and back system

- The seat and back system contains of a huge variety of different upholstered custom-made seat and back units.



### Adjustable head supports

- The different types of head supports are adjustable in angle, height and depth.



### High variety of options

- Other options such as belts, therapy trays or grip ring knobs are available. In addition, individual spoke protectors can be chosen for further personalisation.

#### Frame colors glossy (RAL-colors)



#### Frame colors matt (RAL-colors)



[ottobock.com.au](http://ottobock.com.au)

