



Feature	kidevo mini	kidevo mini.t	kidevo prime	kidevo prime.t
Max. user weight	50 kg	50 kg	75 kg	75 kg
Seat width range	16 – 32 cm	16 – 32 cm	26 – 50 cm	26 – 50 cm
Seat depth range	16 - 35 cm	16 - 35 cm	20 - 50 cm	20 - 50 cm
Tilt	-	-5° to +40°	-	-10° to +50°
Drive wheel sizes	20, 22 and 24°	20, 22 and 24°	22, 24 and 26°	22, 24 and 26°
Crash tested (acc. ISO 7176-19)	Yes	Yes	Yes	Yes
Weight	7,4 kg	9 kg	17,4 kg	19,4 kg